


FTIR bench top spectrometer Interspec 200-X



Specifications:

Wavelength range, IR	7000 to 400 cm^{-1}
Wavelength range, NIR	14000-5800 cm^{-1} or 9000 -3850 cm^{-1}
Resolution, standard	1 cm^{-1}
Resolution, option	0.5 cm^{-1}
Beam at sample	10 mm dia.
Operating system	Windows based
Data acquisition system	18 bit, high speed
Interface	USB 2.0
Power	12 VDC, 30 W
Sample compartment	W20xD26xH16 cm
Dimensions	W59xD39xH19 cm
Weight	24 kg

Interspectrum OU - Toravere, 61602 Tartumaa, Estonia
Tel.: +372 7383008, Fax: +372 7383065
Email: info@interspectrum.ee
www.interspectrum.ee