

## Portable FTIR spectrometer Interspec 301-X



### Specifications:

Wavelength range, KBr	7000 to 400 $\text{cm}^{-1}$
Wavelength range, ZnSe	5000 to 600 $\text{cm}^{-1}$
Resolution, standard	1 $\text{cm}^{-1}$
Resolution, option	0.5 $\text{cm}^{-1}$
Sample compartment	Open architecture
Accessories	Ordinary FTIR accessories
Data acquisition system	18 bit, high speed
Operating system	Windows based
Interface	USB 2.0
Power	12 VDC, 30 W
Dimensions	W49xD39xH20 cm
Weight	18 kg

**Interspectrum OU** - Toravere, 61602 Tartumaa, Estonia  
Tel.: +372 7383008, Fax: +372 7383065  
Email: [info@interspectrum.ee](mailto:info@interspectrum.ee),  
[www.interspectrum.ee](http://www.interspectrum.ee)