

Portable FTIR spectrometer Interspec 300-X



Specifications:

Wavelength range	7000 to 400 cm^{-1} (without ATR accessory)
Wavelength range	6000 to 600 cm^{-1} (with ZnSe ATR accessory)
Resolution, standard	1 cm^{-1}
Resolution, option	0.5 cm^{-1}
Sample compartment	Integrated
Accessories	ATR, reflection, transmission
Data acquisition system	24 bit, high speed
Operating system	Windows based
Interface	USB 2.0
Power	12 VDC, 30 W
Dimensions	W49xD39xH20 cm
Weight	18 kg

Interspectrum OU - Tornii 3, Toravere, 61602 Tartumaa, Estonia
Tel.: +372 7383008, Fax: +372 7383065
Email: info@interspectrum.ee,
www.interspectrum.ee